

IMACTIS[®] CT-Navigation[™] NaviKit[™]

NaviKit is a single-use kit dedicated for percutaneous interventions performed with IMACTIS CT-Navigation. It is used with a non-sterile needle holder for planning.

Needle Holders

The IMACTIS Needle Holder is required to perform a procedure with CT-Navigation. It enables dynamic navigation and virtual trajectory projection in CT images.

There are two types of Needle Holders associated with an IMACTIS CT-Navigation procedure, one is dedicated for sterile conditions, and one is dedicated for planning.

Specific Geometry

- Instrument designed to fit the EM receiver geometry
- Finger placeholder for easy manipulation of the needle holder
- Needle holder axis compatible with linear instruments from 25G to 11G

Blue Non-Sterile Needle Holders

- Provided with NaviKit purchase¹
- Dedicated for planning







Transparent Sterile Needle Holders

- Included inside of NaviKit
- Dedicated for sterile condition
- Single use



Single-Use Components of NaviKit

	Sterile Transparent Adhesive Film	Sterile Needle Holder	Sterile Cover	GripLock
				
Application	For fiducial fixation on patient skin	Single use needle holder	Protection of the EM Sensor	Sterile Cable Securements
Material	Polyurethane with acrylate adhesive	Polycarbonate, Spring: Nickel-cobalt alloy	Latex Free Polymer	N/A

Catalog #	Description
E85101LE	IMACTIS 1 Box of 5 NaviKits + 5 Non-Sterile Needle Holders
E85101LG	IMACTIS 1 Box of 25 NaviKits + 25 Non-Sterile Needle Holders

How to purchase the Navikit

In the U.S.

Call: 866-281-7545

Fax: 877-279-6990

Email: gehccaccessorysales@ge.com

Shop online

www.gehealthcare.com/accessories

Learn more

Contact your local GE HealthCare representative



Navikit is intended for use with CT-Navigation solution. For information regarding the system contact your local GE HealthCare representative.

