

# SIGNA™Works 26.0 Upgrade Order Guide

Welcome to the SIGNA™Works 26.0 software upgrade guide. Please use the steps and information below to configure the appropriate upgrade from your current equipment to the latest performance of SIGNA™Works 26.0 software version.

**If you are a 32-channel system user, select from the Bundled Catalogs or follow the Configuration Steps below. If you are a 16-channel system user, you must follow the Configuration Steps.**

## Bundled Catalogs

If you want the latest technology available, you can also select a pre-bundled upgrade from the list below which will provide all items available.

Catalog	Description	Catalogs Included
<b>1.5T 32CH Systems</b>		
<b>S7526CN</b>	1.5T SIGNA™Works 26.0 Software, Host, and ICN Upgrade	S7526CY, M7005HT
<b>3.0T 32CH Systems</b>		
<b>S7526CP</b>	3.0T SIGNA™Works 26.0 Software, Host, and ICN Upgrade	S7526CY, M7005HT, M7006DA
<b>BrainWave PA and Hardware Users</b>		
<b>S7526CR</b>	BrainWave Upgrade	M7006BF, M7006BG, M7006AK

## Configuration Steps

<b>STEP 1</b>		
Current Equipment	Select	Upgrade
You currently have 8400 host PC or older	➔	SIGNA™Works 26.0 Software and GOC Upgrade ( <b>S7526CY</b> )
You currently have HP z400 or z420 host	➔	SIGNA™Works 26.0 Software and Host Upgrade ( <b>S7526CZ</b> )
You currently have Dell T5180 host	➔	SIGNA™Works 26.0 Software Upgrade ( <b>M7006SW</b> )
<b>STEP 2</b>		
Current Equipment	Select	Upgrade
32CH	You currently have <b>2x ICN Gen2</b>	➔ ICN Upgrade Gen6 32ch Performance ( <b>M7005HT</b> )
	You currently have <b>ICN Gen5, Gen4 or Gen3</b>	➔ ICN Upgrade Gen6 32ch Performance ( <b>M7005HT</b> ) <i>Optional*</i>
16CH	You currently have <b>ICN Gen2</b>	➔ ICN Upgrade Gen6 16ch Performance ( <b>M7005HR</b> )
	You currently have <b>ICN Gen5, Gen4 or Gen3</b>	➔ ICN Upgrade Gen6 16ch Performance ( <b>M7005HR</b> ) <i>Optional*</i>

\* Recommended for enhanced performance. Some advanced features will require specific Generational hardware. ICN considerations:

- **FOCUS** requires Gen4 or higher ICN.
- **DISCO** requires Gen5 or higher 32CH ICN.
- **HyperWorks** and **ViosWorks** require Gen5 or higher ICN. **ViosWorks** also requires a 32-channel system.



STEP 3		
Current Equipment	Select	Upgrade
You currently have the Discovery™ MR750 or MR750w system	➔	3.0T Cable Upgrade Kit <b>(M7006DA)</b>
You currently have the Optima™ MR450 or MR450w system	➔	Nothing Required
STEP 4		
Current Equipment	Select	Upgrade
If you are new to SIGNA™Works 26.0 software, you will need to select a training package:	➔	Apps Training – Base Package (3 Days) <b>(W0138MR)</b>
		Apps Training – Adv Package (4 Days) <b>(W0139MR)</b>
		Apps Training – Adv Package w/Cardiac (6 Days) <b>(W0140MR)</b>
If you are an existing SIGNA™Works 26.0 software customer, training is optional	➔	Training is optional
STEP 5		
Current Equipment	Select	Upgrade
You currently have BrainWave PA	➔	Skull Stripping Upgrade <b>(M7006BF)</b>
In addition to BrainWave PA, if you currently have BrainWave Hardware (M3335MJ)	➔	BrainWave On Host Upgrade <b>(M7006BG)</b> <b>AND</b> BrainWave Hardware Button Response Box <b>(M7006AK)</b>
If you do not have BrainWave PA or hardware	➔	Nothing Required
Order materials complete. Proceed to cart.		

Notes:
1. If you are a CADStream user and not on CADStream 5.2.2, you will need to order additional catalogs. Contact your sales representative.
2. Certain applications like Silent, Silent MRA, Silent PROPELLER and ViosWorks require 32-channel system.

**Need further assistance? Call 1-800-437-1171.**